# **DOMUS AUTO**

## ***Instruction Manual***

**Model No.: 9P-077520**

*PLEASE READ ALL INSTRUCTIONS BEFORE USE.*

# ***USER'S MANUAL***

English

User's manual.....P. 1

Español

Manual de Instrucciones .....P.18

Português

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Français

Instructions d'utilisation.....P.52

German

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Nederlands

Gebruikershandleiding .....P.86

Italiano

Manuale di istruzioni .....P.103

## IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USE

### **DANGER - To reduce the risk of electrical shock:**

1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquids.
5. Do not reach for a product that has fallen into water. Unplug immediately.

### **NOTES, CAUTIONS AND WARNINGS:**

NOTE - Calls attention to useful information, tips, reminders.

CAUTION - Calls attention to correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property.

WARNING - Calls attention to a potential danger that requires correct operating procedures or practices in order to prevent personal injury.

### **WARNING - To reduce the risk of burns, electrical shock, fire, or injury:**

1. This product should never be left unattended when plugged in.
2. Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
3. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.
4. Use this product only for its intended use as described in this manual. Do not use other mattress not recommended by the manufacturer.
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or exposed to water. Return the product to a service center for examination and repair.
6. Keep the cord away from heated surfaces.
7. Never block the air openings of this product. Never place the product on a soft surface, such as a bed or couch, where the air openings may become blocked. Keep the air openings free of lint, hair, and other similar particles.
8. Never drop or insert any object into any opening or hose on this product.
9. Do not modify this equipment without authorization of the manufacturer.
10. Mattress covers have passed skin sensitization and skin irritation test. However, If you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
11. Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.

### **CAUTION –**

1. If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

**SYMBOL Definitions**

Authorized representative in the European community.

Catalog, reorder or reference number

Manufacturer

Complies with standards protecting against electric shock for type BF equipment.

Attention, you should read the accompanying information carefully!

Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

Consult operating instructions for use

**IP21** 

Protected against solid foreign objects of 12,5 mm and greater; Protection against vertically falling water drops

Class II

Temperature limitation/temperature range

Dry clean, Any Solvent Except Trichloroethylene

Do Not Iron

Tumble Dry, Normal, Low Heat

Do Not Tumble Dry

Do Not Bleach

Do Not Dry Clean

Machine wash, regular / normal, 95 degrees C (203 degrees F)

Machine wash, regular / normal, 60 degrees C (140 degrees F)

# 1. INTRODUCTION

This manual should be used for initial set up of the system and for reference purposes.

## 1.1 General Information

This system is a high quality, affordable mattress system, suitable for the treatment and prevention of pressure ulcers. This product has been tested and successfully approved for the following standards:

EN 60601-1  
EN 60601-1-2  
EN 55011 Class B  
IEC61000-3-2  
IEC 61000-3-3

## 1.2 Intended Use

This product is intended:


- to help and reduce the incidence of pressure ulcers while optimizing patient comfort.
- for long term home care of patients suffering from pressure ulcers.
- for pain management as prescribed by a physician.


**NOTE:** This equipment is not suitable for use in the presence of flammable anesthetic mixture with air or with pure oxygen or nitrous oxide.

## 2. PRODUCT DESCRIPTION


### 2.1 Pump and Mattress System

### 2.2 Pump Unit


Front  
1. Quick Connectors  
2. Front Panel

Rear  
1. Power Switch  
2. Power Cord  
3. Hanger  
4. Air Filter

## 2.3 Front panel

### 1. Pressure Adjust

The function can automatically adjust the mattress pressure based on the patient's weight. To remind user the process is on-going, the function LED lights will be periodically flashed from left to right until the Ready LED turn on permanently (the process is completed). There are three conditions will trigger this function to ensure the appropriate pressure setting for the patient.

- A) Press the Auto-Adjust button for 2 seconds, system will automatically re-adjust the pressure setting.
- B) The system will automatically start Auto-Adjust immediately after the initial mattress inflation is completed.
- C) While the system detects a significant change in patient's weight on the mattress (e.g. patient Ingress) for a period of time, the system will automatically run the function again.

While the process is completed, system will return to previous setting automatically. To disable this function, press Auto-Adjust button again.

**NOTE:** During the Auto-Adjust operation, it is normal that system will go through a series of inflation and deflation.

### 2. TUNING

The Tuning function is to fine-tune the system auto-adjusted preset pressure. Three Tuning levels are available. "0" is system auto-adjusted preset pressure. "+" is one level higher than preset pressure, and "-" is one level lower than preset pressure.


**NOTE:** The Tuning function is only available in alternate and static mode.

### 3. THERAPY

The Therapy button allows you to choose one of the therapy modes. There are four selectable therapy modes available. Continually press Therapy button to select a desired mode:

#### A. **MaxFirm**

The surface will go into Max firm automatically when the power switch is turned on. This insures the pump is able to reach its maximum operating pressure. Once the max pressure level is reached, the pump will automatically switch into alternating mode. User can also use this function during patient transfers or nursing procedures for better support. This mode will last for 20 minutes and then it will return to previous setting automatically. To disable this function, press Maxfirm button again.

#### B. **Alternate**

The mattress will alternate every 10 minutes under this therapy mode.

#### C. **Static**

Non-alternating mode, all of the air cells is equally inflated.

#### D. **Seat Inflate**

The seat inflate features additional supports to the patient during upright position.

### 4. Alarm Mute

Press alarm mute button to temporary suspend the low-pressure alarm (LED light and buzzer). Should the situation not resolved within 5 mins, the alarm shall resume to notify the patient.

#### A. **PFA (Power Failure Alarm)**

During power failure situation, the Power Failure LED light will lit on with buzzer. By pressing the mute button to disable both buzzer and LED.

#### B. **Low Pressure Indicator**

When low pressure LED lights up, the pressure inside of air mattress is below normal. Please refer to troubleshooting.

#### C. **Tech. Support**

This feature will light during mechanical failure situation. User can notify the technician for repair.

**5. Panel Lock** Sec

If the control panel remains untouched for 5 minutes, a green LED will light and the control panel will lock all settings. This will prevent the settings from being accidentally changed during normal operation. To unlock, simply press the Panel lock button for three seconds.

### 3. INSTALLATION


Unpack the box and check the package contents for completeness.

**Package Content List**

- Mattress unit x 1 (may not be included if only pump unit is purchased)
- Pump Unit x 1
- User Manual x 1

Inspect equipment for damage which may have occurred during shipment. If there is damage, please contact your dealer immediately.


#### 3.1 Pump & Mattress Installation

1. Place the mattress or pad on top of the bed frame. Please note for the foot end. If it is an overlay system (mattress height less than 8”), please make sure there is a foam or a current mattress underneath.

 **WARNING:** The overlay mattress must be applied on the underlying mattress.

2. Hang the pump onto bed rail (foot-end), and adjust hangers to best upright position of the pump, or place the pump on a flat surface.

3. Connect air hose connectors from air mattress to the pump unit. When a "click" sound is felt or heard, the connection is completed and secured.

 **NOTE:** Check and ensure the air hoses are not kinked or tucked under mattress.

4. Plug the power cord into electrical outlet.

 **NOTE:**

1. Make sure the pump unit is suitable for the local power voltage. Do not position the equipment so that it is difficult to operate the disconnecting device
2. The plug is also used to disconnect the device.

 **CAUTION:** The pump can only be applied to the mattress recommended by the manufacturer. Do not use it for any other purpose. (Applied part: air mattress)

5. Then turn the main power switch to ON position.

 **NOTE:** You can unplug the unit to power-off the device.

 **CAUTION:** The pump supplied with the mattress can only be used for mattresses recommended by the manufacturer. Do not use it for any other purpose. (applied part: air mattress)

 **NOTE:** (For models WITHOUT the low-air-loss function) During power outages, you can cover the Quick Connector with the transport cap to maintain air pressure inside the air cells.

 **NOTE:** After installation, make sure any excess cord or tubing is removed from possible foot traffic to avoid accidental tripping. All EQUIPMENT should be positioned to always allow unhindered patient access by physicians and caregivers.


## 4. OPERATIONS

 **NOTE:** Always read the operating instructions before use.

### 4.1 General operation


# Pumping

1. Switch on the main power switch found on the side of the pump.
2. Once the pump begins delivering air into the mattress it will take approximately 30 minutes to fully inflate the mattress.
3. Every time when mattress is first setup for use, It will be forced to execute Maxfirm for the quickest inflation. The low-pressure indicator (yellow LED) will light up when the mattress is not fully inflated at initial inflation. If the appropriate pressure is reached, the low-pressure indicator (yellow LED) will go off.
4. When the initial inflation (Maxfirm process) is completed, the system will automatically enter the Auto-Adjust process.
5. According to the weight and height of the patient, the system will automatically adjust an appropriate pressure setting to ensure the patient's comfort, and then enter to alternating mode.

### 4.2 Emergency CPR Operations

CPR must be performed on a firm surface. Therefore, if an emergency CPR situation occurs with the patient on the mattress, the mattress must be quickly deflated. This is done by quickly pulling on the CPR tag located at the head of the mattress on the patient's right-side. The quick connector found from the pump unit can also be disconnected to hasten deflation.

### 4.3 Low Pressure Function

When an abnormal low pressure is occurred, the Low Pressure indicator (yellow LED) will light up. Check that connections are correctly made and that they are correctly installed as per installation instructions.

**NOTE** - If the pressure level is consistently low, check for any leakage (tubes or connecting hoses). If necessary, replace any damaged tubes or hoses. Or contact local qualified dealer for repair.

If the pump unit is equipped with alarm buzzer, the alarm will buzz and the Low Pressure indicator light up when air pressure is below normal. To mute the alarm, simply press the  button on the panel. However, the low pressure indicator remains until the low-pressure fault condition is solved.

### 4.4 Static Mode

Press the  button to suspend the alternating mode of mattress. By pressing this button again, the mattress will be back to alternating therapy mode.

### 4.5 Alarm Mute

When low-pressure situation occurred, both the LED light and buzzer will sound off to warn the patient. By pressing the button, it will temporary mute the buzzer so the care giver may check for possible air leaks. Should the situation not resolved within 3 minutes; the alarm will resume to beep.

## 5. CLEANING

It is important to follow the cleaning procedures before first-time use with any patient; this is to avoid the possibility of transmission of pathogens between patients and caregivers.

Wipe down the pump unit with a damp cloth pre-soaked with a mild detergent. Avoid contact with dust and proximity to dusty areas. Make sure that any cleaning agents you use will not harm or corrode the plastic casing on the pump unit.

 **CAUTION:** Do not immerse or soak the pump unit in liquids.

Wipe-down the mattress unit with a damp cloth pre-soaked with warm water containing a mild detergent. Avoid dust and proximity to dusty areas. The cover may also be cleaned using sodium hypochlorite diluted in water. All components should be air dried thoroughly before use.

 **CAUTION:** Do not use phenolic based products for cleaning.

 **CAUTION:** Do not dry the mattress in direct sunlight.

The carrying bag (if available) should be turned inside-out and completely wiped-down using disinfectant solutions. Allow it to air-dry thoroughly and then turn it back and wipe down the outside of the bag with disinfectant solutions.

Cover Material: Nylon/PU	    
Cover Material: Silver+ Nylon/PU	
Cover Material: Stretch	    

## 6. MATTRESS STORAGE

1. Lay the mattress on a flat surface, upside-down.
2. Roll-up the mattress starting from the head until you reach the foot of the mattress.
3. Use the straps at the end of the mattress and stretch them around the mattress to prevent it from unraveling.


 **Note:** Do not fold, crease or stack the mattresses.

## 7. MAINTENANCE

### 7.1 General

1. Check the main power cord and plug for abrasions or excess wear.
2. Check the mattress cover for signs of wear or damage.
3. Disconnect the air tube from the mattress. Then check the airflow coming from the two air outlets on the pump. They should be alternately delivering air when the pump is set at “alternate” mode.
4. Check the air hoses for kinks or breaks. For replacement, please contact your local distributor.

### 7.2 AIR FILTER REPLACEMENT

1. Open the air filter cap located at the back of pump.
2. The filter is reusable and can be washed gently with a mild detergent and water. Dry the filter before use.
3. Check and replace air filter regularly if environment is dirty.

## 8. EXPECTED SERVICE LIFE:

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

## 9. TROUBLE SHOOTING

Problem	Solution
Power is not ON	<ul style="list-style-type: none"> <li>• Check if the plug is connected to mains.</li> </ul>
Alarm is on	<ul style="list-style-type: none"> <li>• Check if the power is suddenly shut down.</li> <li>• Check if the CPR is sealed.</li> <li>• Check if the connection between air tube connector to pump unit is tightly secured.</li> <li>• Check if all tubing connections along mattress are secured.</li> </ul>
Patient is bottoming out	<ul style="list-style-type: none"> <li>• Pressure setting might be inadequate for the patient, adjust comfort range 1 to 2 levels higher and wait for a few minutes for best comfort.</li> </ul>
Mattress form is loose	<ul style="list-style-type: none"> <li>• Check if all the snap buttons or straps of mattress are all securely fastened.</li> <li>• Check if the mattress is fixed to the bed frame by straps.</li> </ul>
No air produced from some air outlets of the air tube connector	<ul style="list-style-type: none"> <li>• This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time.</li> </ul>

## 10. TECHNICAL DESCRIPTION

### SPECIFICATIONS:

Item		Specification
Power Supply (Note: See rating label on the product)		AC 220-240V 50 Hz, 0.08A (for 230V system)
Fuse Rating		T1AL, 250V
Cycle time		Fixed
Dimension (L x W x H)		29 x 18.5 x 12.6 (cm) / 11.4" x 7.3" x 5.0"
Weight		2.3 Kg or 5.1 lb
Environment	Temperature	Operation: 10° C to 40° C (50° F to 104° F) Storage: -15° C to 50° C (5° F to 122° F) Shipping: -15° C to 70° C (5° F to 158° F)
	Humidity	Operation: 10% to 90% non-condensing Storage: 10% to 90% non-condensing Shipping: 10 % to 90% non-condensing
	Atmospheric Pressure	Operation: 700 hPa to 1013.25 hPa Storage: 50 – 106 kPa Shipping: 50 – 106 kPa
Classification		Class II, Type BF, IP21 Applied Part: Air Mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)
Mattress		Specification
Model	5"+ 3" Replacement	8" Replacement
Dimension (L x W x H)	200 x 90 x 12.5 cm / 78.7" x 35.4" x 5"	200 x 90 x 20.3 cm / 78.7" x 35.4" x 8"
Weight	5.6 kg or 12.5 lb	6.9 kg 15.2 lb
Max. Support Weight	180 Kg	200 Kg
Pressure Range	6 comfort levels setting.	

#### NOTE:

1. Consult the distributor or EU representative for further technical documents.
2. These specifications are also applicable for other regions operating with the same power supply.
3. Mattress dimensions and weight is measured without foam;
4. The manufacturer reserves the right to modify the specifications without notice.

## APPENDIX A: EMC INFORMATION

### Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group 1	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

### Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC61000-4-4	±2kV for power supply line ±1kV for input/out line	±2kV for power supply line ±1kV for input/out line	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines	<5 % U <sub>T</sub> (>95 % dip in U <sub>T</sub> )for 0,5 cycle 40 % U <sub>T</sub> (60 % dip in U <sub>T</sub> )for 5 cycles	<5 % U <sub>T</sub> (>95 % dip in U <sub>T</sub> ) for 0,5 cycle 40 % U <sub>T</sub> (60 % dip in U <sub>T</sub> ) for 5 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires

IEC61000-4-11	70 % U <sub>T</sub> (30 % dip in U <sub>T</sub> )for 25 cycles <5 % U <sub>T</sub> (>95 % dip in U <sub>T</sub> )for 5 sec	70 % U <sub>T</sub> (30 % dip in U <sub>T</sub> ) for 25 cycles <5 % U <sub>T</sub> (>95 % dip in U <sub>T</sub> ) for 5 sec	continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U <sub>T</sub> is the a.c. mains voltage prior to the application of the test level			

**Guidance and Manufacturer’s Declaration - Electromagnetic Immunity:**

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance	Electromagnetic Environment-Guidance
Conducted RF IEC 61000-4-6	3Vrms150 kHz to 80 MHz outside ISM bands <sup>a</sup>	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance  <math>d = 1.2\sqrt{P}</math> 150kHz to 80MHz  <math>d = 1.2\sqrt{P}</math> 150kHz to 80MHz  <math>d = 2.3\sqrt{P}</math> 80 MHz to 2.5G MHz</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).<sup>b</sup></p> <p>Field strengths from fixed RF transmitters, as determined by an</p>

		<p>electromagnetic site survey<sup>c</sup>, should be less than the compliance level in each frequency range<sup>d</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	--	--

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a/ The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

b/ The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.

c/ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

d/ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and this device:**

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Print-2014/All rights reserved