Pro-care Auto Pro-care Auto Bariatric Pro-care Auto Pediatric

User's Manual

CONTENTS

ENGLISH

USER MANUAL 1

ESPAÑOL

MANUAL DE INSTRUCCIONES 19

PORTUGUÊS

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FRANÇAIS

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DEUTSCH

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ITALIANO

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DUTCH

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MODEL NO.: 9P-087520

PLEASE READ THE MANUAL BEFORE USE

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USING

DANGER - To reduce the risk of electrocution:

1. Always unplug this product immediately after using.
2. Do not use while bathing.
3. Do not place or store this product where it can fall or be pulled into a tub or sink.
4. Do not place in or drop into water or other liquid.
5. Do not reach for a product that has fallen into water. Unplug immediately.

WARNING - To reduce the risk of burns, electrocution, fire, or injury to persons:

1. Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
2. The product may be used for patients with spinal injury, but consulting a physician before use is suggested. However, it should not be used for patients with unstable spinal fractures.
3. Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.
4. Use this product only for its intended use as described in this manual. Do not use other mattress not recommended by the manufacturer.
5. Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your supplier or Apex Medical Corp. for examination and repair.
6. Keep the cord away from heated surfaces.
7. Never block any air openings of this product or place it on soft surfaces, such as a bed or couch, where openings may be blocked. Keep the air opening free of lint, hair and other similar particles.
8. Never drop or insert any object into any opening or hose.
9. Do not modify this equipment without authorization of the manufacturer.
10. Mattress covers have passed skin sensitization and skin irritation test. However, if you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
11. Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.
12. Third conductor in the POWER SUPPLY CORD is only a functional earth.

CAUTION –

If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

NOTE, CAUTION AND WARNING STATEMENTS:

NOTE - Indicates some tips.

CAUTION - Indicates correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property

WARNING - Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

SYMBOLS

	Authorized representative in the European community.
	Manufacturer
	Complies with standards protecting against electric shock for type BF equipment.
	Consult operating instructions for use
	Attention, should read the instructions.
IP21 	Protected against solid foreign objects of 12,5 mm and greater; Protection against vertically falling water drops
	Class II
	Temperature Limitation
	Dry clean, Any Solvent Except Trichloroethylene
	Do Not Iron
	Tumble Dry, Normal, Low Heat
	Do Not Tumble Dry
	Do Not Bleach
	Do Not Dry Clean
	Machine wash, regular / normal, 95 degrees C (203 degrees F)
	Machine wash, regular / normal, 60 degrees C (140 degrees F)
	Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product.

1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

1.1 General Information

The system is a high quality and affordable mattress system suitable for treatment and prevention of pressure ulcers.

The system has been tested and successfully approved to the following standards:

IEC/EN 60601-1
IEC/EN 60601-1-2
IEC/EN 61000-3-2 Class A
IEC/EN 61000-3-3
CISPR 11 Group 1, Class B

EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.


1.2 Intended Use

This product is intended to help and reduce the incidence of pressure ulcers while optimizing patient comfort. It also provides following purposes:

- to help and reduce the incidence of pressure ulcers while optimizing patient comfort.
- for long term care of patients suffering from pressure ulcers.
- for pain management as prescribed by a physician.

The product can only be operated by personnel who are qualified to perform general nursing procedures and have received adequate training in knowledge of prevention and treatment of pressure ulcers.

 **NOTE:** The equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.


 **WARNING:** Do not turn or flip patients who have bone fracture or spinal cord injury.

2. Product Description

Unpack the box to check for any damage which may have occurred during shipment. If there is any damage, please contact the point of purchase immediately.

2.1 Pump & Mattress System


1. C.P.R.
2. Pump Unit
3. Mattress system
4. Quick Connector

2.2 Pump Unit


Front

1. Quick Connector
2. Front Panel


Side

1. Power Cord
2. Fuses


Rear

1. Mounting brackets
2. Air filter

2.3 Front Panel

1. Power button 
2. Pressure Adjustment 

The function can automatically adjust the mattress pressure based on the patient's weight. To remind user the process is on-going, the function LED lights will periodically flash and be constantly on after the process is completed. There are three conditions will trigger this function to ensure the appropriate pressure setting for the patient.

- A) When in therapy mode, press the TUNING button for 3 seconds, system will automatically start the readjustment of the pressure setting.
- B) The system will automatically start Auto-Adjust immediately after the initial mattress inflation is completed.

C) While the system detects a significant change in patient's weight on the mattress (e.g. patient Ingress) for a period of time, the system will automatically run the function again.

NOTE: Every time the mattress is initialized (inflated), it will automatically go to "Max Firm" mode to hasten inflation. Once the system is ready to use, the system will go to alternate mode automatically.


NOTE: You can also check if the pressure is suitable for the patient by sliding one hand beneath the air cells at the level of the patient's buttocks. Always leave at least 1" inch space between patient and the static cell to prevent bottoming out.


NOTE: During the Auto-Adjust operation, it is normal that system will go through a series of inflation and deflation.

3. TUNING

The Tuning function is to fine-tune the system auto-adjusted preset pressure. Three Tuning levels are available. "0" is system auto-adjusted preset pressure. "+" is one level higher than preset pressure, and "-" is one level lower than preset pressure.

 **NOTE:** The Tuning function is available in Alternating, Alternating Low Pressure and Continuous Low Pressure mode.

 **NOTE:** Always perform a hand-check after pressure setting to prevent the patient from bottoming out.

4. Cycle Time

There are four selectable cycle times can be found on the control panel. By pressing the button, user can select one of four cycle times based upon patient comfort and desired outcome.

5. Modes

Each button allows you to choose one of the modes. There are four selectable modes available. Press the button to select the desired mode:

5.1 Alternate Mode


Alternating

Continuously and sequentially inflates and deflates air cells to achieve periodic pressure relief and avoid prolonged pressure on any single point of the patient's body.

Alternating
Low
Pressure

Continuously and sequentially inflates and deflates air cells to achieve periodic pressure relief and avoid prolonged pressure on any single point of the patient's body. The air cells are inflated at lower pressure and deflated at higher pressure when compared to the respective comfort level in alternating mode.

 **NOTE:** The air cells in the head section of Pro-care Auto Pediatric system are also alternating.

5.2 Continuous Low Pressure Mode

Constantly non-alternating mode. The system will redistribute body mass over a greater surface area at a constant low pressure. All of the air cells are equally inflated at lower pressure when compared to the respective comfort level in alternating mode.

5.3 Max Firm

The surface will go into Max Firm automatically when the power switch is turned on. This ensures the pump is able to reach its maximum operating pressure. Once the max pressure level is reached, the pump will automatically switch into Auto-Adjust mode. User can also use this function during patient transfers or nursing procedures for better support. This mode will last for 20 minutes and then it will return to previous setting automatically. To disable this function, press Max Firm button again.

 **NOTE:** Max Firm and preset function mode LED indicator will be on simultaneously during initial mattress inflation until the mattress reaches its maximum pressure.

6. Panel Lock

If the control panel remains untouched for 5 minutes, a green LED will light and the control panel will lock all settings. This will prevent the settings from being accidentally changed during normal operation. To unlock, simply press the panel lock button for **three** seconds.

7. Alarm Mute

When the Alarm Mute button is pressed, the buzzer will be temporarily suspended. If the situation is not resolved within 3 minutes, the alarm will be activated again.

7.1 PFA (Power Failure Alarm)

During a power failure situation, the Power Failure LED light and buzzer will be activated. To disable acoustic alarm, press the mute button.

7.2 Low Pressure Indicator

When low pressure LED lights up, the pressure inside of air mattress is below normal. Please refer to troubleshooting.


7.3 High Pressure Indicator

When high pressure lights up, the pressure inside of air mattress is higher than normal. Please refer to troubleshooting.

7.4 Tech. Support

The Service indicator light will be activated when an issue occurs within the system. User can notify the authorized technician for repair.


3. Installation

- ⚠ NOTE:** Unpack the box to inspect for any damage which may have occurred during shipment. If there is any damage, please contact your dealer immediately.
- Place the mattress on top of the bed frame. There are securing straps on the base of the mattress. Secure the mattress firmly by fixing the straps to the bed frame; ensure that moving sections of the bed frame are still free to move.
 - ⚠ WARNING:** The 5" (12.7cm) overlay mattress must be applied on the underlying mattress.
 - Hang the pump onto bed rail (foot-end), and adjust hangers to best upright position of the pump, or place the pump on a flat surface.
 - Connect air hose connectors from air mattress to the pump unit. A "click" sound will be heard when connection is completed and secured.
 - ⚠ NOTE:** Check and ensure the air hoses are not kinked or tucked under mattress.
 - Plug the power cord into electrical outlet.
 - ⚠ CAUTION:** Make sure the pump unit is suitable for the local power voltage and frequency.
 - ⚠ CAUTION:** Position the device in a place where it is easy to disconnect.
 - Then turn the main power switch to ON position.
 - ⚠ CAUTION:** The pump should only be used with the mattress recommended by the manufacturer. Do not use it for any other purpose.
 - The cable management loops are provided along the left side of the mattress to cover the power cable. Please make sure the cable has no signs of damage. After installation, the extra length of the power cord, if any, should be neatly managed to avoid causing a trip and should be clear of moving bed mechanisms or other possible entrapment area.
- The EQUIPMENT should be firmly placed at position where users/doctors can access to easily.

4. Operation

⚠ NOTE: Always read the operating instructions before use.

4.1 General Operation

1. Switch on the main power switch found on the side of the pump.
2. Push the power button to turn on and run.
3. The system will automatically go into Max Firm mode for the quickest inflation. Once the pump begins delivering air into the mattress it will take approximately 15-30 minutes to fully inflate the mattress.

⚠ NOTE: Please do not place patient on a mattress until it is fully inflated.

4. During the initial inflation (Max Firm process), user may preset the system by pressing the Mode button and select a desired mode, while both selected mode and Max Firm LED indicators will be simultaneously on.
5. When the initial inflation (Max Firm process) is completed, the system will automatically enter the Auto-Adjust process. During the Auto-Adjust process, the preset function LED indicator will be off.
6. According to patient weight differences, the system will automatically adjust an appropriate pressure setting to ensure the patient's comfort, and then go into the alternating mode or the previously set mode.

⚠ NOTE: The last five cells from the foot-end have heel-relief connectors. Carers can easily disconnect the connector and heels can be suspended in air to completely eliminate pressure exerted on the heels.


4.2 Pressure FINE-TUNE

Depending on the desired preference, user can adjust mattress pressure slightly by pressing Tuning button .


⚠ NOTE: Check to see if the suitable pressure is selected by sliding one hand between the deflated air cells and the patient to feel patient's buttocks. Users should be able to feel the minimum contact.

4.3 Emergency CPR Operations

- ⚠ NOTE:** When there is an emergency to perform CPR on the patient, quickly turn the CPR knob to release air from mattress. The CPR knob is located at the head-end, right-hand-side of the mattress. The quick connector on the pump unit can be disconnected for even faster deflation.
- ⚠ NOTE:** After CPR, re-activate auto adjustment to ensure proper pressure setting. Please make sure the CPR valve and quick connector are correctly connected, and always performs hand-check after Auto-Adjust or Tuning process.

4.4 Transportation

When there is a power outage or a need of transporting the patient, the quick connector can be capped as shown to slow down the air leakage. Depending on the setting and mattress height, the mattress can remain inflated for a period of time.

- ⚠ WARNING:** If there is a bottom-out issue, please remove the patient from deflating mattress to another suitable surface.

5. Cleaning

It is important to follow the cleaning procedures to avoid cross contamination. Be sure to clean the surface in a dry and dust free environment.

A) Pump Unit

Wipe down the pump unit with a damp cloth pre-soaked with a mild detergent. Avoid contact with dust and proximity to dusty areas. Make sure that any cleaning agents you use will not harm or corrode the plastic casing on the pump unit. If your doctor or medical facilities have other special cleaning instructions, please follow the professional instructions.

 **CAUTION:** Do not immerse or soak pump unit.

 **WARNING:** Do not remove the housing of the pump to avoid electrical shock. All disassembly or repair should be done by professional technicians.

 **CAUTION:** The pump does not need oil lubrication; please do not disassemble the system.

B) Mattress Cover

Cover Material : Stretch	    
Cover Material: Nylon PU	    

Wipe-down the mattress unit with a damp cloth pre-soaked with warm water containing a mild detergent, or chlorine bleach followed by an approved intermediate level disinfectant. Also the mattress top cover can be completely removed for laundry with water temperature indicated on the laundry label; however, it is recommended that the user still check with local policy to determine the time/ temperature ratio required to achieve thermal disinfection. The cover may also be cleaned using 6% of 1,000 ppm Sodium Hypochlorite diluted in water. After cleaning, please avoid dust and proximity to dusty areas and all parts should be air dried thoroughly before use.

 **CAUTION:** Do not use phenolic based products for cleaning.

 **CAUTION:** After cleaning, dry the mattress without direct exposure to sunlight.

6. Storage

1. Disconnect the quick connector
2. Roll from the foot end towards the head end with CPR valve open, and make sure the tubing is not kinked
3. Foot-end strap can then be stretched around the rolled mattress to prevent unrolling. Fasten the strap buckle to secure the packed mattress.


 **NOTE:** Do not kink, crease or stack the mattresses and do not store the system in direct sunlight, high temperature or moist place.

7. Maintenance

7.1 General


1. Check main power cord and plug for any abrasion or excessive wear.
2. Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are connected together correctly.
3. Check airflow from the quick connector. The airflow should alternate between each connector every half-cycle time if it is in alternating mode.
4. Check the air hoses for kink or breaks. For replacement, please contact our local provider.

7.2 Fuse Replacement

1. Disconnect the plug from mains power when a blown fuse is suspected.
2. Use a proper tool to remove the cover of the fuse holder.
3. Insert a new fuse of the correct rating in, and replace the fuse holder back to correct position. The fuse should be rated as T1.6AL / 250V type and VDE approved.

7.3 Air Filter Replacement

1. Replace the air filter located at the back of pump.
2. The filter is reusable and can be washed gently with a mild detergent and water. Air dry the filter before use.
3. Check and replace air filter regularly.

8. Expected Service Life

The products are intended to offer safe and reliable operation when used or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

9. Trouble Shooting

Q.1 Power is not ON

- Check if the power plug is connected to the power outlet.
- Check for a blown fuse.

Q.2 Low Pressure Alarm is on

- Check if the quick connector is tightly secured.
- Check if the CPR knob is closed.
- Check if all tubing connections along mattress are secured.
- Check if the air hoses are kinked or broken.

Q.3 High Pressure Alarm is on

- Check if the tubing between mattress and pump set is kinked, twisted or bended.

Q.4 Power Failure Alarm is on

- Check if the power is suddenly shut down.
- Check if the power cord is connected properly.

Q.5 Patient is bottoming out

- Please check if the function mode setting is correct.
- Press Tuning button for three seconds to readjust.
- Adjust pressure by using Tuning button + and wait for a few minutes. Perform a hand check again to see if the problems were solved.

Q.6 Mattress form is loose

- Check if all the snap buttons or straps of mattress are all securely fastened.
- Check if the mattress is fixed to the bed frame by straps.

Q.7 No air produced from some air outlets of the air tube connector

- This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time.

If the above information does not solve your problems, please contact the point of purchase or authorized distributor directly. There might require a technician to take care of the problem.

10. Technical Specifications

Control Unit		Specification			
Power Supply (Note: See rating label on the product)		AC 220-240V 50 Hz, 0.15A			
Fuse Rating		T1.6AL, 250V			
Cycle time		Four selectable (10, 15, 20, 25 mins)			
Dimension (L x W x H)		34.0 x 13.5 x 20.5 cm / 13.4" x 5.3" x 8.1"			
Weight		3.8 Kg / 8.37 lb			
Environment	Temperature	Operation: 10°C to 40°C (50°F to 104°F) Storage: -15°C to 50°C (5°F to 122°F) Shipping: -15°C to 70°C (5°F to 158°F)			
	Humidity	Operation: 10% to 90% non-condensing Storage: 10% to 90% non-condensing Shipping: 10 % to 90% non-condensing			
	Atmospheric Pressure	Operation: 700 hPa to 1060 hPa Storage: 700 hPa to 1060 hPa Shipping: 700 hPa to 1060 hPa			
Classification		Class II, Type BF, IP21 Applied Part: Air Mattress Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)			
Mattress					
Model	Dimensions			Weight (Approx).	Max. Weight Capacity
	Length	Width	Height		
5" Overlay	200 cm 78.7"	80/85/90 cm 31.5"/33.5"/35.4"	13 cm 5"	5.9 kg 13 lb	200 kg 440 lb
5"+3" Replacement			20 cm 5"+3"		
8" Replacement	200 cm 78.7"	75/80/85/90/100 cm 29.5"/31.5"/33.5"/35.4"/39.4"	20 cm 8"	10 kg 22 lb	250 kg 550 lb
8" Replacement	210 cm 82.7"	85/90/100 cm 33.5"/35.4"/39.4"	20 cm 8"	10 kg 22 lb	250 kg 550 lb
8" Bariatric Replacement	200 cm 78.7"	107 /122 cm 42"/48"	20 cm 8"	13 kg 28.5 lb	450 kg 990 lb
Pediatric	140/150 cm 55"/59"	65/70 cm 25.5"/27.5"	18 cm 7"	9.5 kg 20.9 lb	95 kg 209 lb

NOTE :

1. Consult the distributor or EU representative for further technical documents.
2. The specifications are also suitable for other areas operating with the same power supply.
3. Mattress dimension and weight is measured without foam base.
4. The manufacturer reserves the right to modify the specifications without prior notice.

Appendix A: EMC Information

Guidance and Manufacturer's Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment-Guidance
RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / Flicker emissions IEC61000-3-3	Complies	

Warning:

1. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.
2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

Basic EMC standard	Immunity Test Levels		Compliance Levels	Electromagnetic Environment-Guidance
	Professional healthcare facility environment	HOME HEALTHCARE ENVIRONMENT		
Electrostatic Discharge (ESD) IEC61000-4-2	±8 kV contact ± 2, 4, 8, 15 kV air		±8 kV contact ± 2, 4, 8, 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC61000-4-4	±2 kV, 100 kHz for power supply line ±1 kV, 100 kHz for input/output line		±2 kV, 100 kHz for power supply line ±1 kV, 100 kHz for input/output line	Mains power quality should be that of atypical commercial or hospital environment
Surge IEC61000-4-5	± 0.5, 1 kV (line to line) ± 0.5, 1, 2 kV (line to earth)		± 0.5, 1 kV (line to line)	Mains power quality should be that of atypical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	Voltage Dips: i) 100% reduction for 0.5 period, ii) 100% reduction for 1 period, iii) 30% reduction for 25/30 period, Voltage Interruptions: 100% reduction for 250/300 period		230 Vac input	Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	30 A/m	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms 0,15 MHz – 80 MHz 6 Vrms in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to the frequency of the transmitter.

<p>Radiated RF EM Fields IEC61000-4-3</p>	<p>3 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz</p> <p>380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation , (1kHz sine)</p>	<p>10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz</p> <p>380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation , (1kHz sine)</p>	<p>10V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz</p> <p>380 - 5800 MHz, 9-28 V/m, Pulse modulation or FM ±5 kHz deviation, (1kHz sine)</p>	<p>Recommended separation distance</p> $d = \sqrt{P}$ <p>150kHz to 80MHz</p> $d = 0.6\sqrt{P}$ <p>80MHz to 800MHz</p> $d = 1.2\sqrt{P}$ <p>800 MHz to 2.7GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency ranged.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1: U_i is the a.c. mains voltage prior to the application of the test level</p> <p>NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people</p>				
<p>a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.</p> <p>b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.</p>				

Recommended separation distances between portable and mobile RF communications equipment and this device:

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \sqrt{P}$	80 MHz to 800 MHz $d = 0.6\sqrt{P}$	800 MHz to 2.7 GHz $d = 1.2\sqrt{P}$
0.01	0.1	0.06	0.12
0.1	0.31	0.19	0.38
1	1	0.6	1.2
10	3.1	1.9	3.8
100	10	6	12

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Print-2021/All rights reserved