# Optima Auto Link

User Manual

Be well, Live well  
**Wellell**

# CONTENTS

## **ENGLISH**

|             |   |
|-------------|---|
| USER MANUAL | 1 |
|-------------|---|

## **ESPAÑOL**

|                         |    |
|-------------------------|----|
| MANUAL DE INSTRUCCIONES | 27 |
|-------------------------|----|

## **PORTUGUÊS**

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## **FRANÇAIS**

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## **DEUTSCH**

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## **DUTCH**

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## **ITALIANO**

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**Model No.:** 9P-086520

PLEASE READ ALL INSTRUCTIONS BEFORE USE

## **IMPORTANT SAFEGUARDS**

### **READ ALL INSTRUCTIONS BEFORE USING**

#### **DANGER - To reduce the risk of electrocution:**

- (1). Always unplug this product immediately after using.
- (2). Do not use while bathing.
- (3). Do not place or immerse the unit into water or other liquids.
- (4). Do not reach for a product that has fallen into water. Unplug immediately.

#### **WARNING - To reduce the risk of burns, electrocution, fire, or injury to persons:**

- (1). Evaluate patients for entrapment risk according to protocol and monitor patients appropriately.
- (2). The product may be used with caution for patients with spinal injury, but suggested to consult with physician before use. However, it should not be used for patients with unstable spinal fractures.
- (3). Close supervision is necessary when this product is used on or near children. Electrical burns or choking accident may result from a child swallowing a small part detached from the device.
- (4). Use this product only for its intended use as described in this manual. Do not use other mattress not recommended by the manufacturer.
- (5). Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your supplier or Wellex Inc. for examination and repair.
- (6). Keep the cord away from heated surfaces.
- (7). Never block any air openings of this product or place it on soft surfaces, such as a bed or couch, where openings may be blocked. Keep the air opening free of lint, hair, and other similar particles.
- (8). Never drop or insert any object into any opening or hose.
- (9). No modification of this equipment is allowed unless it is conducted by the authorized technician.
- (10). Mattress covers have passed skin sensitization and skin irritation test. However, If you suspect that you may have had or are having an allergic reaction, please consult a physician immediately.
- (11). Do not leave long lengths of tubing around the top of your bed. It could lead to strangulation.
- (12). Keep the pump away from flammable liquids or gases.
- (13). Do not use the extension cord for a long time.
- (14). When disposing the unit because of damage or failure, you must comply with the regulations of your particular region or country for disposing of such equipment.
- (15). Medical electrical equipment needs to special precautions regarding EMC and needs to be installed according to EMC information. Contact your service provider or manufacturers regarding EMC installation information.
- (16). If there is a possibility of electro-magnetic interference with mobile phones, please increase the distance (3.3m) between devices or turn off the mobile phone.

**CAUTION – Reducing the risk of damage to the device.**

- (1). Mobile RF communications equipment can affect medical electrical equipment.
- (2). Make sure the filter area at the back of the device is not blocked.
- (3). Do not place the device directly onto flammable materials.
- (4). The equipment must not be covered or positioned in such a way that the operation or performance of the equipment is adversely affected.
- (5). Be careful not to place the device where it can be taken or touched by children and pets.
- (6). Store the product in a dry place. Keep out of direct sunlight. Use of or store the product only in the following conditions:
  - Operation: 10°C to 40°C (50°F to 104°F); 10% to 90% non-condensing
  - Storage: -15°C to 50°C (5°F to 122°F); 10% to 90% non-condensing
- (7). Connection of OUR PRODUCT to 4G LTE NETWORK that includes other equipment could result in previously unidentified RISKS to PATIENTS, OPERATORS or third parties. The USER should identify, analyze, evaluate and control these RISKS. Subsequent changes to the 4G LTE NETWORK could introduce new RISKS and require additional analysis. Changes to the 4G LTE NETWORK include the following:
  - changes in the 4G LTE NETWORK configuration;
  - connection of additional items to the 4G LTE NETWORK;
  - disconnecting items from the 4G LTE NETWORK;
  - update of equipment connected to the 4G LTE NETWORK; and
  - upgrade of equipment connected to the 4G LTE NETWORK.

**NOTE, CAUTION AND WARNING STATEMENTS:**

**NOTE:** Indicate some tips.

**CAUTION:** Indicate correct operating or maintenance procedures in order to prevent damage to or destruction of the equipment or other property

**WARNING:** Calls attention to a potential danger that requires correct procedures or practices in order to prevent personal injury.

**SYMBOLS**

|                                                                                               |                                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                   |                                                                                                 |
|-----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
|              | Authorized representative in the European community.                                                                                                                                                                                                                                                                                                                                                              |  | Temperature Limitation                                                                          |
|              | Manufacturer                                                                                                                                                                                                                                                                                                                                                                                                      |  | Dry clean, Any Solvent Except Trichloroethylene                                                 |
|              | Complies with standards protecting against electric shock for type BF equipment.                                                                                                                                                                                                                                                                                                                                  |  | Do Not Iron                                                                                     |
|              | Consult operating instructions for use                                                                                                                                                                                                                                                                                                                                                                            |  | Tumble Dry, Normal, Low Heat                                                                    |
|              | Attention, should read the instructions.                                                                                                                                                                                                                                                                                                                                                                          |  | Do Not Tumble Dry                                                                               |
|              | Attention, should read the instructions.                                                                                                                                                                                                                                                                                                                                                                          |  | Do Not Bleach                                                                                   |
|              | Do Not Dry Clean                                                                                                                                                                                                                                                                                                                                                                                                  |  | Machine wash, regular / normal, 95 degrees C (203 degrees F)                                    |
|              | Class II                                                                                                                                                                                                                                                                                                                                                                                                          |  | EC directive medical device label (Compliance to MDD 93/42/EEC; 4 digits number is NB's number) |
|              | FCC logo                                                                                                                                                                                                                                                                                                                                                                                                          |                                                                                   |                                                                                                 |
| <b>IP21</b>  | Protected against solid foreign objects of 12.5 mm and greater; Protection against vertically falling water drops                                                                                                                                                                                                                                                                                                 |                                                                                   |                                                                                                 |
|            | Attention – Observe proper Disposal of Electrical & Electronic Equipment (WEEE): This product should be handed over to an appropriate collection point for the recycling of electrical and electronic equipment. For more detailed information about the recycling of this product, please contact your local city office, household waste disposal service or the retail store where you purchased this product. |                                                                                   |                                                                                                 |

# 1. Introduction

This manual should be used for initial set up of the system and for reference purposes.

## 1.1. General Information

The system is a high quality mattress system suitable for treatment and prevention of pressure injury. This is a micro-processor controlled device. It features the illuminated LCD display.

The system has been tested and successfully approved to the following standards:

IEC/EN 60601-1  
IEC/EN 60601-1-2  
IEC/EN 61000-3-2 Class A  
IEC/EN 61000-3-3  
CISPR 11 Group 1, Class B

## EMC Warning Statement

This equipment has been tested and found to comply with the limits for medical devices to the IEC/EN 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

## 1.2. Intended Use

This product is intended to help and reduce the incidence of pressure injuries while optimizing patient's comfort. It can support patient up to 250 kg (550 lb). And it can also provide the following purposes:

- For long term home care of patients suffering from pressure injury.
- For pain management as prescribed by a physician.

The product can only be operated by personnel who are qualified to perform general nursing procedures and have received adequate training in knowledge of prevention and treatment of pressure injury.

 **NOTE:** Equipment is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.

## 2. Product Description

Unpack the box to check for any damage which may have occurred during shipment. If there are damages, please contact your dealer immediately.

### 2.1. Contents

This system comes with the followings items.

Unpack this box by removing the device and its accessory and checking for any damage, which may have occurred during shipping. If there are damages, please contact your service provider immediately.

- |                    |                 |
|--------------------|-----------------|
| (1). Pump          | (4).Air filter  |
| (2). Ac power cord | (5).User manual |
| (3). Mattress      |                 |

 The SIM card is pre-installed and bound to your device. Contact your service provider if it requires replacement.


 NOTE: The content of standard package might vary or change without advanced notification.

 CAUTION: Patient should not connect the device to a personal computer for data downloading. This may cause device failure.

## 2.2. Pump and Mattress System


- ❶ CPR
- ❷ Pump unit
- ❸ Mattress system
- ❹ Quick connector
- ❺ Position sensor

## 2.3. Pump Unit


Front

- ❶ Quick connector
- ❷ Position sensor
- ❸ Front panel


Rear

- ❶ Mounting brackets
- ❷ USB cover
- ❸ SIM card slot
- ❹ Air filter
- ❺ Battery cover
- ❻ Power socket
- ❼ Fuse
- ❽ Power switch

## 2.4. Front Panel

### (1). Power Key

Power Key is to turn on and turn off pump.

### (2). Panel Lock key

User interface will be locked automatically 5 minutes after last operation. To unlock interface, press and hold Panel Lock key for 1 second.

**NOTE:** Always unlock user interface before turning off the pump.

### (3). Max Firm key

This key is to make sure mattress can be reached maximal operating pressure. This mode will last for 20 minutes and then return to previous setting automatically. This function can be used during patient transfers or nursing procedures for better support. To disable this function, please press Max firm key again.

### (4). LCD & Knob

Knob is to choose and confirm setting displayed on LCD, rotate the knob and press once to select the setting mode desired in the highlight.

### (5). Alarms and alarm mute key

When there is an alert situation happened, the LCD screen will display a corresponding error code with an audible buzzer in an interval of a few seconds. The buzzer can be temporarily suspended by pressing the Alarm Mute key. Alarm buzzer will signal repeatedly every 3 minute until the alert situation is resolved.

Below presents 2 types of alarm and please refer to the attached error code table for details.

### (6). Signal icon

This appears ONLY on the IoT-technology-enabled device. The icon shows signal strength. Make sure to place the device at where the signal coverage is sound.


**Note:** The IoT (Internet of Things) function is ON automatically when the system is activated.

## 3. Installation

### 3.1. Pump & Mattress

**⚠ NOTE:** Unpack the box to inspect for any damage which may have occurred during shipment. If there are any damages, please contact your dealer immediately.

- (1). Place the mattress on top of bed frame and foot print toward foot end. Secure the mattress firmly by fixing the straps to bed frame to ensure mattress will be free from sliding. Securing straps are located on the base of the mattress

**⚠ NOTE:** Make sure that ALL the Relief knobs and emergency CPR knob are closed and locked in position.

**⚠ CAUTION:** Placing another layer of sheet or bed linen between the patient and the mattress may reduce the benefits provided by the mattress and should be avoided at all time.

- (2). Hang the pump onto bed panel at the foot-end and adjust hangers to hang the pump to best upright position, or place the pump on a flat surface.
- (3). Connect the air hose connector of air mattress and the electrical cord of position sensor to the pump unit. Make sure the air hose connection is completed and secured when a "click" sound is felt or heard and position sensor assemble the alignment with the arrow pointing up.

**⚠ NOTE:** Check and ensure the air hose is not kinked or tucked under mattress.


**⚠ NOTE:** Auto recognition mattress of Optima Auto.


- (4). Plug the power cord into electrical outlet.
-  **NOTE:** Make sure the pump unit is suitable for the local power voltage.
-  **NOTE:** The plug is also served to stop the device.
-  **NOTE:** Pump might still be in operation after the power cord disconnection if optional battery is built-in.
- (5). Turn the power switch to ON position.
-  **CAUTION:** The pump can only be applied to the mattress recommended by the manufacturer. Do not use it for any other purpose. (Applied part: air mattress) Do not position the equipment so that it is difficult to operate the disconnecting device.
-  **CAUTION:** Allow the mattress to inflate before positioning the patient on the mattress.
- (6). The cable management tape is provided along the side of the mattress to manage power cord. Please make sure the wire has no signs of damage. After installation, extra length of the power cord, if any, should be neatly managed to prevent stumble accident and should be clear of moving bed mechanisms or other possible entrapment areas.
-  **NOTE:** The EQUIPMENT should be firmly placed at position where users can easily access to.

3.2. Position Sensor- installation & removal

Position sensor is key to auto seat inflation. It allows the mattress automatically to inflate when the head section has been raised to prevent bottoming out.

To initiate this function, you should connect the sensor cord to the unit.

| To                                                                                                                                                                                                                                     | Steps                                                                                                                                                                                                                                                                                                                                                                                    | Illustration                                                                                                                                                                      |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Attach                                                                                                                                                                                                                                 | <div><div>(1). Align the arrows on the cord and the unit.</div><div>(2). Insert the cord accordingly as illustrated.</div></div>                                                                                                                                                                                                                                                         |                                                                                                 |
| Detach                                                                                                                                                                                                                                 | <div><div><div>(1). Put your fingers at the end of the cord as illustrated.</div><div>(2). Pull backward to detach it from the unit.</div></div><div> <b>Note:</b> The grey area is the protection mechanism that prevents misuse and enables a proper detachment by sliding it backward.</div></div> | <div></div> |
| <div> <b>Warning:</b> Do not immerse in water.</div> <div></div> |                                                                                                                                                                                                                                                                                                                                                                                          |                                                                                                                                                                                   |

# 4. Operation

 NOTE: Always read the instructions before use.

## 4.1. General Operation

Switch the power switch on the side of the pump and press power key to turn it on.


### 4.1.1. Pre-set program

- (1). Choose “TICK” and press the knob to enter the pre-set program, and then the pump will run itself automatically.


Perform pre-set program:  
Mode → Alternating  
Cycle time → 10 mins.

- (2). Enter setting to set your prefer Mode and Cycle time.

- (3). Once the pump identifies the size of the mattress via the quick connector, it starts to inflate the mattress. The initial inflation takes about 30 minutes.

- (4). After initial inflation, the system sets off to detect the patient's weight in order to provide corresponding pressures for better mattress pressure relieving performance.


The 5 blue bars displayed on the screen represents the status of the automatic adjustment detection. When completes, the 5 bars will stop running and become .

 **NOTE:** Place the patient on the mattress **ONLY** after the initial inflation is completed, better when the auto adjustment is about to start.


 **NOTE:** Once the patient is on the mattress, the pump will automatically sense and adjust the pressure in the cells to support the patient.

 **NOTE:** Press the knob for at least 3 seconds to quickly initiate the auto adjustment process.


- (5). The current therapy mode is shown the LCD display. Turn the knob clockwise or counter-clockwise to change settings. Press the knob when your prefer selection is highlighted to enter or confirm the setting.


#### A. Mode

Alternating


Continuous Low Pressure


Back to previous page

 **NOTE:** Review the patient's monitoring and reposition program before you change the operation mode.

## B. Cycle time

## C. Tuning

**NOTE:** The pressure in the air cells can be manually adjusted for patient comfort by entering the tuning mode in the lower right corner of the LCD screen.

**NOTE:** Slide one hand underneath the air cell where the patient's buttock locates to ensure that the given air pressure is appropriate. Always leave a space for about 2.54 cm (1 inch) between the patient and the cell to prevent bottoming-out.

## 4.2. Internet of Things Service (IoT)

The Optima Auto Link is an IoT-technology-enabled device. This IoT device streams the pump or mattresses' real-time conditions via a cloud-based platform named KeepWell.

For the IoT device, the LCD displays a signal icon that shows the current strength.

For more detail, please refer to the manual of the cloud-based platform:


"KeepWell—streamline operation, logistics, and customer service for your business."

**NOTE:** Use this device in places where 4G LTE NETWORK coverage is good to ensure optimal IoT services.

**NOTE:** Consult with the service provider or the manufacturer for the cellular communication service provider.


## 4.3. Emergency CPR Operation

CPR knob is located at the patient's left-hand side of the mattress near the head area. When there is a need to perform CPR on the patient, quickly turn the CPR knob to release air from mattress. The quick connector on the pump unit can be disconnected for even faster deflation.


#### 4.4. Relief Function

The 5 knobs are on the patient's left-hand side at the foot end of the mattress. Turning the knob to deflate the cell enables the caregivers to float the patient's heels and take away the pressure from bony prominences.

Closed

Open

For permanent off-loading/pressure relief:

- Select **no more than two** cells directly under the area that desired to relieve.
- Rotate the knob in a counterclockwise direction to the 12 o'clock position to deflate the cell.

**⚠ NOTE:** This one to two cells can be left permanently deflated.


**⚠ NOTE:** Deflating more than two adjacent cells may affect the support of the patient during the normal alternating cycle.

#### 4.5. Evacuation Handles

**⚠ NOTE:** When an emergency occurs, the designated handles can be used for caregivers to quickly and safely evacuate the patient from the bed in the facility.

#### 4.6. Transportation

**⚠ NOTE:** When there is a need of patient transportation, pump with the power cord can be disconnected from the mattress during transport. For optimum comfort, WELLELL recommends caregivers to change to Continuous Low Pressure mode and wait for 5 minutes before disconnecting. Cap on the quick connector as shown on the left to keep the air inside the mattress. The mattress remains inflated at least 24 hours depending on the setting. Optional battery model can support at least 6 hours of treatment mode.

#### 4.7. Battery/filter replacement

The device requires battery/filter replacement. Your service provider shall provide the replacement service.

|                           |                           |
|---------------------------|---------------------------|
| The filter: every 6 month | The battery: every 2 year |
|---------------------------|---------------------------|

## 4.8. Built-in batteries

For model with built-in batteries, it provides non-stop 6-hour therapy when in power outage. You can refer to top-right side of LCD display for battery status. Please read information below for battery operation.

- (1). Use AC power whenever possible. Connect the pump to AC power outlet and turn on the power switch to fully charge the battery for emergencies.
- (2). The warranty period is 6 months for the built-in batteries.
- (3). Connect the pump to AC power outlet immediately when the low-battery graphic shows up.
- (4). The battery may not be fully charged upon receipt. Connect the pump to AC power outlet for at least 10 hours.

 **NOTE:** Pump can not complete initial inflation when battery capacity is under 25%.

## 4.9. Cellular communication

The IoT device uses the 4G LTE NETWORK that transmits the device's data to your service provider wirelessly.

- (1). The 4G LTE NETWORK is activated when the system is turned on by default.
- (2). The icon displayed at the top left of your screen indicates signal strength. Use this device in places where 4G LTE NETWORK coverage is good to ensure optimal IoT services.
- (3). The data is automatically transmitted to the service provider. In order to make sure that your data is transferred, leave your device connected to the power outlet at all times.

 **NOTE:** Data might not be transmitted if you use it outside of the country or region of purchase.

 **NOTE:** Devices with cellular communication might not be available in all regions.

 **NOTE:** Please be aware that within the wireless network, the availability and quality of the network may be affected by terrain, buildings, and the weather.

 **NOTE:** Consult with the service provider or the manufacturer for the cellular communication service provider.

# 5. Cleaning

It is important to follow the cleaning procedures to avoid cross contamination. Be sure to clean the surface in a dry and dust free environment. Wipe down the pump unit with a damp cloth pre-soaked with a mild detergent. Avoid contact with dust and proximity to dusty areas. Make sure that any cleaning agents you use will not harm or corrode the plastic casing on the pump unit. If your doctor or medical facilities have other special cleaning instruction, please follow the professional instruction.

-  **CAUTION:** Do not immerse or soak pump unit/position sensor. 
-  **WARNING:** Do not remove the casing of the pump to avoid electrical shock. All disassembly or repair should be done by authorized technicians.
-  **CAUTION:** The pump does not need oil lubrication; please do not dissemble the system.

|                                                                            |                                                                                                                                                                                                                                                                                                                                                                                                                           |
|----------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Cover Material: <ul style="list-style-type: none"><li>• Stretch</li></ul>  |      |
| Cover Material: <ul style="list-style-type: none"><li>• Nylon/PU</li></ul> |      |

Wipe-down the mattress unit with a damp cloth pre-soaked with warm water containing a mild detergent, or chlorine bleach followed by an approved intermediate level disinfectant. Also the mattress top cover can be completely removed for laundry with water temperature indicated on the laundry label; however, it is recommended that the user still check with local policy to determine the time/temperature ratio required achieving thermal disinfection. The cover may also be cleaned using sodium hypochlorite diluted in water. After cleaning, please avoid dust and proximity to dusty areas and all parts should be air dried thoroughly before use.

-  **CAUTION:** Do not use phenolic based products for cleaning.
-  **CAUTION:** Keep mattresses out of direct exposure to sunlight, and make sure to dry the mattress after cleaning.
-  **CAUTION:** Make sure to execute the cleaning procedure every time before a new patient lying on the air mattress.

# 6. Storage

- (1). To store the mattress, lay the mattress out flat and upsides down.
- (2). Roll from the head end towards the foot end with CPR valve open.
- (3). Foot-end strap can then be stretched around the rolled mattress to prevent unrolling.
- (4). After the product is stored for 3 months, it needs to be recharged again if optional battery is built-in.
-  **NOTE:** Do not kink, crease or stack the mattresses and do not store the system in direct sunlight, high temperature or moist areas.


## 7. Maintenance

### 7.1. General

- (1). Check main power cord and plug if there is abrasions or excessive wear.
- (2). Check mattress cover for signs of wear or damage. Ensure mattress cover and tubes are stubbed together correctly.
- (3). Check airflow from the quick connector. The airflow should alternate between two ventilation tubes in the middle and the other two at the bottom of quick connector every half-cycle time if it's in alternating mode.
- (4). Check the air hoses if there is kink or break. For replacement, please contact local representative or place of purchase.


### 7.2. Fuse Replacement

- (1). Disconnect the power plug when a blown fuse is suspected.
- (2). Remove the cover of the fuse holder by means of a flat screwdriver.
- (3). Insert a new fuse of the correct rating in, and replace the cover of the fuse holder back. The fuse should be rated as T1AL/250V type and VDE approved.

### 7.3. Air Filter Replacement

- (1). Open the air filter cap located at the back of pump.
- (2). Replace with a clean filter and put the cover lid back. The filter is reusable and can be washed gently with a mild detergent and water. Dry the filter before use.
- (3). Check and replace air filter regularly if environment is dirty. Suggested to do it **every six month**.

## 8. Expected Service Life

The products are intended to offer safe and reliable operation when use or installed according to the instructions provided by Apex Medical. Apex Medical recommends that the system be inspected and serviced by authorized technicians if there are any signs of wear or concerns with device function and indication on products. Otherwise, service and inspection of the devices generally should not be required.

## 9. Trouble Shooting

### 9.1. User tackle-able issues

Issues that users might encounter and solve by themselves. See below for the description, corresponding displays and solutions.

#### (1). Low Pressure Alarm

Alarm will be indicated when the pump detects low pressure condition within the mattress.



 NOTE: In alternating mode, the system first inflates and then has the cells hold at a certain level to support the patient while monitoring the cells at the same time. The situation continues until the low pressure alarm is lifted.


 NOTE: The pump seals the air inside the mattress, which temporarily keeps the patients from bottoming out during low pressure condition. Patient should be moved if the power failure issue cannot be solved soon.

#### (2). Low battery capacity

 NOTE: When low battery symbol is shown on the display, plug the power cord into the wall outlet immediately.

 CAUTION: Pump system will stop running; the LCD displays the alarm, along with the buzzer images when the battery capacity is blow 10%. To restore the system, plug the power cord into the wall outlet, and press the power key to restart the pump.

| Code No.   | Diagram                                                                             | Description                      | What to Do                                                                                                                                                                                                 |
|------------|-------------------------------------------------------------------------------------|----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 01, 02, 05 |    | Low pressure alarm: body section | (1).Check all the joint parts to see if they are well connected, or if there is leakage, along with kinked or broken parts that need to be addressed.<br><br>(2).Call for service if the problem persists. |
| 03, 04, 06 |   | Low pressure alarm: foot section |                                                                                                                                                                                                            |
| 07         |  | Low pressure alarm: head section |                                                                                                                                                                                                            |
| 19         |  | Position sensor is not connected | Re-start the pump to fix this. Call for service if the problem persists.                                                                                                                                   |

| Code No. | Diagram                                                                            | Description                                                                                                      | What to Do                                                                                                                                                                                    |
|----------|------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 20       |   | Quick connector is not connected                                                                                 | (1). Check if the joint parts are well connected. Ensure that there is no dirty or foreign objects blocking the quick connector.<br>(2). Replace the quick connector if the problem persists. |
| 21       |   | Low battery capacity (< 25 %)                                                                                    | Plug the power cord into the wall outlet, and then re-start the pump.                                                                                                                         |
| 21_1     |   | Low battery capacity (< 10 %)                                                                                    |                                                                                                                                                                                               |
| 21_2     |   | Initial inflation:<br>Low battery capacity (< 10 %)<br>After initial inflation:<br>Low battery capacity (<= 2 %) |                                                                                                                                                                                               |
| 24       |  | Connect to a wrong mattress                                                                                      | Ensure you have attached the pump to the right mattress.                                                                                                                                      |

## 9.2. Technician tackle-able issues

Pump issues that need to be addressed by professional technicians. The LCD displays each error code and an accompanying cogwheel-with-wrench symbol to alert you. If encountering this, **contact service center for further assistance.**

### (1). Power failure alarm

-  NOTE: There will be buzz sound and the alarm mute LED will light up constantly with dark screen when POWER FAILURE occurs.
-  NOTE: The pump can seal the air inside the mattress, which temporarily keeps the patients from bottoming out during power failure. Patients should be moved if the power failure issue cannot be solved soon.

| Code No. | Description                         | What to Do                                                                                                                                                                                                                                                      |
|----------|-------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 00       | Initial inflation overtime          | <ul style="list-style-type: none"> <li>•To initiate, battery capacity must be up to 25%.</li> <li>•Check all the joint parts to see if they are well connected, or if there is leakage, along with kinked or broken parts that need to be addressed.</li> </ul> |
| 08-16    | System malfunction                  | Call customer service for further assistance.                                                                                                                                                                                                                   |
| 17       | Pressure control system malfunction | <ul style="list-style-type: none"> <li>•Re-start the system.</li> <li>•Re-connect the pressure sensor.</li> <li>•Make sure the tube is not kinked or broken.</li> <li>•Call for service if the problem persists.</li> </ul>                                     |
| 18       | System access control malfunction   | Call customer service for further assistance.                                                                                                                                                                                                                   |

## (2). Communication failure

| Code No. | Description                           | What to Do                                                                                                                                                                                                                       |
|----------|---------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 60       | SIM card not found                    | Call customer service for further assistance.                                                                                                                                                                                    |
| 61       | Cellular module communication failure |                                                                                                                                                                                                                                  |
| 62       | Cellular module activation failure    |                                                                                                                                                                                                                                  |
| 63       | Data upload failure                   |                                                                                                                                                                                                                                  |
| 64       | Server communication failure          | (1). Ensure that the 4G LTE NETWORK signal at the top left corner on the display is strong.<br>(2). If not, relocate the system to a place where the cellular signal is sound.<br>(3). Call for service if the problem persists. |
| 65       | Weak LTE signal                       | Call customer service for further assistance.                                                                                                                                                                                    |

### Caution

- Relevant risks need to be analyzed and evaluated by manufacturer. Relevant protective measures should be taken.
- Any changes to the relevant network configuration and network device connection methods may cause unknown risks. If necessary, please contact your service provider or the manufacturer.

### 9.3. Others

**Q1 Power is not ON**

- Check if the plug is connected to mains.
- Check for a blown fuse.

**Q2 Low Pressure Alarm is on**

- Check if the Quick Connector is tightly secured.
- Check if all tubing connections along with mattress are secured.
- Check if the air hoses are kinked or broken.
- Check if the CPR knob is tightly closed.

**Q3 Power Failure Alarm is on**

- Check if the power is suddenly shut down.
- Check if the power cord is connected properly.

**Q4 Patient is bottoming out**

Pressure setting might be inadequate for the patient, increase comfort range 1 level higher and wait for a few minutes for best comfort.

**Q5 Mattress form is loose**

- Check if all the snap buttons or straps of mattress are all securely fastened.
- Check if the mattress is fixed to the bed frame by straps.

**Q6 No air produced from some air outlets of the air tube connector**

This is normal since there is alternating mode. Air outlets take turns to produce air during their cycle time.

**Q7 System malfunction (Code 11 or 12)**

Check if there is only one or two knobs at the foot-end of the mattress remained open.

If the above information does not solve your problems, please contact your local agent directly. They might require a technician to take care of the problem.

# 10. Technical Specifications

| Item                                                 |                      | Specification                                                                                                                                                                                                                                                                                                                                                                                                 |
|------------------------------------------------------|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Power Supply (Note: See rating label on the product) |                      | AC 220-240V, 50/60Hz, 0.6A (with battery)                                                                                                                                                                                                                                                                                                                                                                     |
| Fuse Rating                                          |                      | T1AL, 250V                                                                                                                                                                                                                                                                                                                                                                                                    |
| Cycle time                                           |                      | 10, 15, 20, 25 minutes                                                                                                                                                                                                                                                                                                                                                                                        |
| Dimension (L x W x H)                                |                      | 34.1 x 16.5 x 26.0 cm                                                                                                                                                                                                                                                                                                                                                                                         |
| Weight                                               |                      | 5 Kg /11 lbs (with battery)                                                                                                                                                                                                                                                                                                                                                                                   |
| Environment                                          | Atmospheric Pressure | 700 hPa to 1060 hPa                                                                                                                                                                                                                                                                                                                                                                                           |
|                                                      | Temperature          | Operation: 10°C to 40°C (50°F to 104°F)<br>Storage: -15°C to 50°C (5°F to 122°F)<br>Shipping: -15°C to 70°C (5°F to 158°F)                                                                                                                                                                                                                                                                                    |
|                                                      | Humidity             | Operation: 10% to 90% non-condensing<br>Storage: 10% to 90% non-condensing<br>Shipping: 10 % to 90% non-condensing                                                                                                                                                                                                                                                                                            |
| Classification                                       |                      | Class II, Type BF, IP21<br>Applied Part: Air Mattress<br>Not suitable for use in the presence of a flammable anesthetic mixture (No AP or APG protection)                                                                                                                                                                                                                                                     |
| Communication                                        |                      | <b>Cellular</b><br>*Twelve Bands FDD-LTE: 700, 800, 850, 900, 1700/2100 (AWS), 1800, 1900, 2100, 2600 MHz (bands 1, 2, 3, 4, 5, 7, 8, 12, 18, 19, 20, 28)<br>* Seven Bands UMTS (WCDMA/FDD): 800, 850, 900, 1700/2100 (AWS), 1800, 1900 and 2100 MHz (bands 1, 2, 4, 5, 8, 9, 19)<br>*Quad Band GSM: 850, 900, 1800 and 1900 MHz<br><br>Bd1, Bd2, Bd3, Bd4, Bd5, Bd7, Bd8, Bd12, Bd17, Bd18, Bd19, Bd20, Bd28 |
|                                                      |                      | RF ID No: FCC ID: QIPPLS62-W<br> <b>Note: The device should be used at a minimum distance of 8" (20 cm) from the body during operation.</b>                                                                                                                                                                                |

| Mattress              | Specification                                            |
|-----------------------|----------------------------------------------------------|
| Model                 | 8" Mattress                                              |
| Dimension (L x W x H) | 200 x 80/85/90 x 20 cm<br>78.7" x 31.5"/33.5"/35.4" x 8" |
| Weight                | 13 Kg / 28.6 lbs                                         |
| Max. Support Weight   | 250 Kg / 550 lbs                                         |
| Relief Function       | 5 knobs(16 <sup>th</sup> to 21 <sup>st</sup> air cells)  |

**NOTE:**

- Consult the distributor or EU representative for further technical documents.
- The specifications are also suitable for other areas operating with same power supply.
- Mattress dimension and weight are measured without foam base.
- The manufacturer reserves the right to modify the specifications without notice.
- Cycle time may vary according to the local power supply.

# Appendix A: EMC Information

## Guidance and Manufacturer’s Declaration- Electromagnetic Emissions:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

| Emissions Test                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | Compliance | Electromagnetic Environment-Guidance                                                                                                                                           |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions<br>CISPR 11                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | Group1     | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment |
| RF emissions<br>CISPR 11                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | Class B    | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network            |
| Harmonic emissions<br>IEC61000-3-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Class A    |                                                                                                                                                                                |
| Voltage fluctuations /<br>Flicker emissions<br>IEC61000-3-3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | Complies   |                                                                                                                                                                                |
| <div> Warning:</div> <div><div>1.</div><div>The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.</div></div> <div><div>2.</div><div>Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.</div></div> <div><div>3.</div><div>Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.</div></div> |            |                                                                                                                                                                                |

### Guidance and Manufacturer's Declaration- Electromagnetic Immunity:

This device is intended for use in the electromagnetic environment specified below. The user of this device should make sure it is used in such an environment.

| Basic EMC standard                                                                                    | Immunity Test Levels                                                                                                                                                                          |                                                                                              | Compliance Levels                                                    | Electromagnetic Environment-Guidance                                                                                                                                                                                                                                          |
|-------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|----------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                                       | Professional healthcare facility environment                                                                                                                                                  | HOME HEALTHCARE ENVIRONMENT                                                                  |                                                                      |                                                                                                                                                                                                                                                                               |
| Electrostatic Discharge (ESD)<br>IEC61000-4-2                                                         | ±8 kV contact<br>± 2, 4, 8, 15 kV air                                                                                                                                                         |                                                                                              | ±8 kV contact<br>± 2, 4, 8, 15 kV air                                | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.                                                                                                                                |
| Electrical fast transient/ burst<br>IEC61000-4-4                                                      | ±2 kV, 100 kHz for power supply line<br>±1 kV, 100 kHz for input/output line                                                                                                                  |                                                                                              | ±2 kV, 100 kHz for power supply line                                 | Mains power quality should be that of atypical commercial or hospital environment                                                                                                                                                                                             |
| Surge<br>IEC61000-4-5                                                                                 | ± 0.5, 1 kV (line to line)<br>± 0.5, 1, 2 kV (line to earth)                                                                                                                                  |                                                                                              | ± 0.5, 1 kV (line to line)                                           | Mains power quality should be that of atypical commercial or hospital environment.                                                                                                                                                                                            |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC61000-4-11 | Voltage Dips:<br>i) 100% reduction for 0.5 period,<br>ii) 100% reduction for 1 period,<br>iii) 30% reduction for 25/30 period,<br>Voltage Interruptions:<br>100% reduction for 250/300 period |                                                                                              | 230 Vac                                                              | Mains power quality should be that of a typical commercial or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60Hz) magnetic field<br>IEC61000-4-8                                              | 30 A/m                                                                                                                                                                                        | 30 A/m                                                                                       | 30 A/m                                                               | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                     |
| Conducted RF<br>IEC 61000-4-6                                                                         | 3 Vrms<br>0,15 MHz – 80 MHz<br>6 Vrms in ISM bands between 0,15 MHz and 80 MHz<br>80 % AM at 1 kHz                                                                                            | 3 Vrms<br>0,15 MHz – 80 MHz<br>6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 | 3 Vrms<br>0,15 MHz – 80 MHz<br>6 Vrms in ISM and amateur radio bands | Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there                                                                                                                                             |

|                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                  |                                                                                                                                                     |                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                           |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                  | MHz<br>80 % AM at 1 kHz                                                                                                                             | between<br>0,15 MHz and<br>80 MHz<br>80 % AM at 1<br>kHz                                                                                                        | commended separation<br>distance calculated from<br>the equation applicable to<br>the frequency of the<br>transmitter.                                                                                                                                                                                                                                                                                                    |
| Radiated RF EM<br>Fields<br>IEC61000-4-3                                                                                                                                                                                                                                                                                                   | 3 V/m<br>80 MHz - 2.7 GHz<br>80 % AM at 1 kHz<br><br>380 - 5800 MHz,<br>9-28 V/m,<br>Pulse modulation or<br>FM ±5 kHz deviation ,<br>(1kHz sine) | 10 V/m<br>80 MHz - 2.7 GHz<br>80 % AM at 1 kHz<br><br>380 - 5800 MHz,<br>9-28 V/m,<br>Pulse modulation<br>or FM ±5 kHz<br>deviation, (1kHz<br>sine) |                                                                                                                                                                 | Recommended separation<br>distance<br>$d=\sqrt{P}$ 150kHz to 80MHz<br>$d=0.6\sqrt{P}$ 80MHz to<br>800MHz<br>$d=1.2\sqrt{P}$ 800 MHz to<br>2.7GHz                                                                                                                                                                                                                                                                          |
|                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                  |                                                                                                                                                     | 10 V/m<br>80 MHz - 2.7<br>GHz<br>80 % AM at 1<br>kHz<br><br>380 - 5800<br>MHz,<br>9-28 V/m,<br>Pulse<br>modulation<br>or FM ±5 kHz<br>deviation,<br>(1kHz sine) | Where P is the maximum<br>output power rating of the<br>transmitter in watts (W)<br>according to the<br>transmitter manufacturer<br>and d is the<br>recommended separation<br>distance in meters (m). <sup>b</sup><br><br>Field strengths from fixed<br>RF transmitters, as<br>determined by an<br>electromagnetic site<br>survey , <sup>a</sup> should be less<br>than the compliance level<br>in each frequency ranged. |
|                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                  |                                                                                                                                                     |                                                                                                                                                                 | Interference may occur in<br>the vicinity of equipment<br>marked with the following<br>symbol:<br><br>                                                                                                                                                                                                                                 |
| NOTE 1: U <sub>T</sub> is the a.c. mains voltage prior to the application of the test level<br>NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.<br>NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by<br>absorption and reflection from structures, objects and people |                                                                                                                                                  |                                                                                                                                                     |                                                                                                                                                                 |                                                                                                                                                                                                                                                                                                                                                                                                                           |

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and this device:**

This device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of this device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and this device as recommended below, according to the maximum output power of the communications equipment

| Rated maximum output power of transmitter<br>W                                                                                                                                                                                                                                                                                          | Separation distance according to frequency of transmitter    m |                                        |                                         |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|----------------------------------------|-----------------------------------------|
|                                                                                                                                                                                                                                                                                                                                         | 150 kHz to 80 MHz<br>$d = \sqrt{P}$                            | 80 MHz to 800 MHz<br>$d = 0.6\sqrt{P}$ | 800 MHz to 2.7 GHz<br>$d = 1.2\sqrt{P}$ |
| 0.01                                                                                                                                                                                                                                                                                                                                    | 0.1                                                            | 0.06                                   | 0.12                                    |
| 0.1                                                                                                                                                                                                                                                                                                                                     | 0.31                                                           | 0.19                                   | 0.38                                    |
| 1                                                                                                                                                                                                                                                                                                                                       | 1                                                              | 0.6                                    | 1.2                                     |
| 10                                                                                                                                                                                                                                                                                                                                      | 3.1                                                            | 1.9                                    | 3.8                                     |
| 100                                                                                                                                                                                                                                                                                                                                     | 10                                                             | 6                                      | 12                                      |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance <i>d</i> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. |                                                                |                                        |                                         |
| Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.                                                                                                                                                                                                                                          |                                                                |                                        |                                         |
| Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.                                                                                                                                                                    |                                                                |                                        |                                         |

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## Guidance and Manufacturer’s Declaration- Electromagnetic Immunity:

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| Basic EMC standard                            | Immunity Test Levels                         |                             |  | Compliance Levels                     | Electromagnetic Environment-Guidance                                        |
|-----------------------------------------------|----------------------------------------------|-----------------------------|--|---------------------------------------|-----------------------------------------------------------------------------|
|                                               | Professional healthcare facility environment | HOME HEALTHCARE ENVIRONMENT |  |                                       |                                                                             |
| Electrostatic Discharge (ESD)<br>IEC61000-4-2 | ±8 kV contact<br>± 2, 4, 8, 15 kV air        |                             |  | ±8 kV contact<br>± 2, 4, 8, 15 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with |

|                                                                                                       |                                                                                                                                                                                               |                                                                                                                      |                                                                                                  |                                                                                                                                                                                                                                                                              |
|-------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                                       |                                                                                                                                                                                               |                                                                                                                      |                                                                                                  | synthetic material, the relative humidity should be at least 30 %.                                                                                                                                                                                                           |
| Electrical fast transient/ burst<br>IEC61000-4-4                                                      | ±2 kV, 100 kHz for power supply line<br>±1 kV, 100 kHz for input/output line                                                                                                                  |                                                                                                                      | ±2 kV, 100 kHz for power supply line                                                             | Mains power quality should be that of atypical commercial or hospital environment                                                                                                                                                                                            |
| Surge<br>IEC61000-4-5                                                                                 | ± 0.5, 1 kV (line to line)<br>± 0.5, 1, 2 kV (line to earth)                                                                                                                                  |                                                                                                                      | ± 0.5, 1 kV (line to line)                                                                       | Mains power quality should be that of atypical commercial or hospital environment.                                                                                                                                                                                           |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC61000-4-11 | Voltage Dips:<br>i) 100% reduction for 0.5 period,<br>ii) 100% reduction for 1 period,<br>iii) 30% reduction for 25/30 period,<br>Voltage Interruptions:<br>100% reduction for 250/300 period |                                                                                                                      | 230 Vac                                                                                          | Mains power quality should be that of a typical commercial or hospital environment. If the use of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60Hz) magnetic field<br>IEC61000-4-8                                              | 30 A/m                                                                                                                                                                                        | 30 A/m                                                                                                               | 30 A/m                                                                                           | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                    |
| Conducted RF<br>IEC 61000-4-6                                                                         | 3 Vrms<br>0,15 MHz – 80 MHz<br>6 Vrms in ISM bands between 0,15 MHz and 80 MHz<br>80 % AM at 1 kHz                                                                                            | 3 Vrms<br>0,15 MHz – 80 MHz<br>6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz<br>80 % AM at 1 kHz | 3 Vrms<br>0,15 MHz – 80 MHz<br>6 Vrms in ISM and amateur radio bands between 0,15 MHz and 80 MHz | Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than there commended separation distance calculated from the equation applicable to                                                                   |

|                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                               |                                                                                                                                            |                                                                                                                                            |                                                                                                                                                                                                                                                                                                                                                                                    |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                               |                                                                                                                                            | 80 % AM at 1 kHz                                                                                                                           | the frequency of the transmitter.                                                                                                                                                                                                                                                                                                                                                  |
| Radiated RF EM Fields<br>IEC61000-4-3                                                                                                                                                                                                                                                                                                   | 3 V/m<br>80 MHz - 2.7 GHz<br>80 % AM at 1 kHz<br><br>380 - 5800 MHz,<br>9-28 V/m,<br>Pulse modulation or FM ±5 kHz deviation ,<br>(1kHz sine) | 10 V/m<br>80 MHz - 2.7 GHz<br>80 % AM at 1 kHz<br><br>380 - 5800 MHz,<br>9-28 V/m,<br>Pulse modulation or FM ±5 kHz deviation, (1kHz sine) | 10 V/m<br>80 MHz - 2.7 GHz<br>80 % AM at 1 kHz<br><br>380 - 5800 MHz,<br>9-28 V/m,<br>Pulse modulation or FM ±5 kHz deviation, (1kHz sine) | Recommended separation distance<br>$d=\sqrt{P}$ 150kHz to 80MHz<br>$d=0.6\sqrt{P}$ 80MHz to 800MHz<br>$d=1.2\sqrt{P}$ 800 MHz to 2.7GHz                                                                                                                                                                                                                                            |
|                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                               |                                                                                                                                            |                                                                                                                                            | Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). <sup>b</sup><br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey , <sup>a</sup> should be less than the compliance level in each frequency ranged. |
|                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                               |                                                                                                                                            |                                                                                                                                            | Interference may occur in the vicinity of equipment marked with the following symbol:<br><br>                                                                                                                                                                                                   |
| NOTE 1: U <sub>r</sub> is the a.c. mains voltage prior to the application of the test level<br>NOTE 2: At 80 MHz and 800 MHz, the higher frequency range applies.<br>NOTE 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people |                                                                                                                                               |                                                                                                                                            |                                                                                                                                            |                                                                                                                                                                                                                                                                                                                                                                                    |

- a)Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

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| Rated maximum output<br>power<br>of transmitter<br>W                                                                                                                                                                                                                                                                               | Separation distance according to frequency of transmitter    m |                                      |                                       |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|--------------------------------------|---------------------------------------|
|                                                                                                                                                                                                                                                                                                                                    | 150 kHz to 80 MHz<br>$d=\sqrt{P}$                              | 80 MHz to 800 MHz<br>$d=0.6\sqrt{P}$ | 800 MHz to 2.7 GHz<br>$d=1.2\sqrt{P}$ |
| 0.01                                                                                                                                                                                                                                                                                                                               | 0.1                                                            | 0.06                                 | 0.12                                  |
| 0.1                                                                                                                                                                                                                                                                                                                                | 0.31                                                           | 0.19                                 | 0.38                                  |
| 1                                                                                                                                                                                                                                                                                                                                  | 1                                                              | 0.6                                  | 1.2                                   |
| 10                                                                                                                                                                                                                                                                                                                                 | 3.1                                                            | 1.9                                  | 3.8                                   |
| 100                                                                                                                                                                                                                                                                                                                                | 10                                                             | 6                                    | 12                                    |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. |                                                                |                                      |                                       |
| Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.                                                                                                                                                                                                                                     |                                                                |                                      |                                       |
| Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.                                                                                                                                                               |                                                                |                                      |                                       |

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Print-2023-03-13 /All rights reserved
**Apex Medical Corp.**