

WIZARD 310/320

Nasal mask/ Full face mask

Flexibility and comfort for better daily sleep

Individuals with sleep-disorder breathing can experience improvements in their sleeping quality as quickly as two weeks of positive airway therapy and significantly reduce daytime sleepiness and its impact on quality of life. Factors such as properly minimizing breathing obstruction incidents, improving mask comfort, and reducing the risk of pressure sores, air leakage around the cushion, claustrophobia, nasal congestion, and additional complications from frequently treatment sessions drastically impacts therapy compliance.^{1,2}

WIZARD 310&320's ergonomic cushion optimizes pressure redistribution, pain-free forehead support, non-magnetic snap-on and -off parts, and tubing safety valve improve ease-of-use and comfort, allowing users to integrate CPAP into their daily routine with confidence quickly.

1. Patel SR, White DP, Malhotra A, Stanchina ML, Ayas NT. Continuous Positive Airway Pressure Therapy for Treating gess in a Diverse Population With Obstructive Sleep Apnea: Results of a Meta-analysis. Arch Intern Med. 2003;163(5):565-571. doi:10.1001/archinte.163.5.565.
2. Kohler M, Smith D, Tippet V, Stradling JR. Predictors of long-term compliance with continuous positive airway pressure. Thorax. 2010 Sep;65(9):829-32. doi: 10.1136/thx.2010.135848. PMID: 20805182.



Ultimate fit

The ergonomic cushion with multiple-pressure-redistribution design provides better facial support to each individual.



Comfortable

Flexible and pain-free frames allow more mobility during leisure activities and sleep, improving the therapy experience.

Simplifying sleep therapy, improving sleep quality and long-term therapy compliance.



Soft Forehead Support

The soft and flexible forehead support structure with less skin contact provides further comfort and retains stability.



Flexible Mask Frame

The adaptable and lightweight frame offers better fitting and great seal.



Breathable and Durable Headgear

Ensures a comfortable and effective OSA treatment.




Safety Valve Structure

Designed to avoid asphyxia in case of CPAP machine power failure (for WiZARD 320 only)



Quick Release Tubing


An easy one step removal design for convenient tubing connection. Offers better flexibility for freedom of movement during sleep.



Snap-on design



Buckle system



Comfortable all night

Soft side frames are gentle on the face and cheekbones, ensuring comfort even during long periods of wear.



Quiet

Exhaling at a 45° angle reduces noise and air leaks, helping better sleep for the user and their partner.



Simple

Parts that are easy to assemble, store, and clean motivate users to stay consistent in their therapy when traveling.

